

G gauge	OD	length				
		10mm	13mm	15mm	16mm	25mm
15G	1.80mm			30°		
17G	1.40mm		30°			
18G	1.20mm		30°			
19G	1.10mm		30°			
20G	0.90mm	45°	45° 30°		45°	
21G	0.80mm		45° 0°			
22G	0.70mm		45° 30° 0°			
23G	0.60mm		45° 30° 0°			
25G	0.50mm		45° 0°			
27G	0.40mm		45° 30° 0°			
30G	0.30mm					0°

* We are able to handle anything from supply of our standard size injection needles on an OEM basis to customer requirements for specified dimensions.

* Manufacture of special sizes of injection needles is accommodated through sample manufacture.