

美容分野 美容針

Aesthetic Field / Aesthetic Needles

Gゲージ	外径	長さ											
		4mm	6mm	8mm	12mm	16mm	25mm	30mm	40mm	50mm	75mm	100mm	125mm
■ 15G	1.80mm						●						
□ 16G	1.60mm						●						
■ 17G	1.40mm							●					
■ 18G	1.20mm						●			●	●	●	●
■ 19G	1.10mm												
■ 20G	0.90mm									●			
■ 21G	0.80mm						●			●			
■ 22G	0.70mm							●					
■ 23G	0.60mm						●			●			
■ 24G	0.55mm						●						
■ 25G	0.45mm						●			●			
■ 26G	0.40mm				●	●		●		●			
■ 27G	0.40mm	●	●	●	●	●	●			●			
■ 30G	0.30mm	●	●	●	●		●		●				
■ 31G	0.25mm	●	●	●	●		●						
■ 32G	0.23mm	●	●	●	●								
■ 33G	0.20mm	●	●	●	○								
■ 34G	0.18mm	○	○	○									

社の標準サイズの注射針のOEM提供からお客様の要求寸法まで対応可能です。

●特注サイズの注射針の製造は、サンプル製造から対応いたします。

● STANDARD SIZES

○ UNDER REQUEST SIZES