

ホスピタル分野 皮下注射針

Hospital Field / Hypodermic Needles

Gゲージ	外径	長さ																
		4mm	6mm	10mm	12mm	16mm	19mm	25mm	30mm	38mm	50mm	60mm	75mm	80mm	90mm	100mm	120mm	125mm
■ 14G	2.10mm							○	○	○	○	○	○	○	○	○	○	○
■ 15G	1.80mm							○	○	○	○	○	○	○	○	○	○	○
□ 16G	1.60mm						○	○	○	○	○	○	○	○	○	○	○	○
■ 17G	1.40mm						○	○	○	○	○	○	○	○	○	○	○	○
■ 18G	1.20mm				●	●	●	●	●	●	●	○	○	○	○	○	○	○
■ 19G	1.10mm				○	●	●	●	●	●	●	○	○	○	○	○	○	○
■ 20G	0.90mm			○	●	●	●	●	●	●	●	○	○	○	○	○	○	○
■ 21G	0.80mm			○	●	●	●	●	●	●	●	○	○	○	○	○	○	○
■ 22G	0.70mm			○	○	○	●	●	●	●	●	○	○	○				
■ 23G	0.60mm			○	○	●	●	●	●	●	●	○	○	○				
■ 24G	0.55mm			○	○	○	○	●	●	●	○	○						
■ 25G	0.45mm			○	○	●	●	●	●	●	●	○						
■ 26G	0.40mm			○	●	●	●	●	●	●								
■ 27G	0.40mm	○	○	○	●	●	●	●	●	●								
■ 28G	0.36mm	○	○	○	●	○	○	○	○									
■ 29G	0.33mm	○	○	○	○	○	○	○	○									
■ 30G	0.30mm	○	○	○	○	○	○											
■ 32G	0.23mm	○	○	○	○													
■ 33G	0.20mm	○	○	○	○													

社の標準サイズの注射針のOEM提供からお客様の要求寸法まで対応可能です。

●特注サイズの注射針の製造は、サンプル製造から対応いたします。

● STANDARD SIZES
○ UNDER REQUEST SIZES