
## **Closed End**



	Ggauge	OD	length
	<b>2</b> 1G	0.80mm	23mm
	<b>2</b> 3G	0.60mm	23mm
	25G	0.50mm	23mm
]	27G	0.40mm	23mm
	30G	0.30mm	23mm

## **Blunt Tip**

	Ggauge	OD	length
0	20G	0.90mm	22mm
	<b>25</b> G	0.50mm	40mm
	<mark>3</mark> 0G	0.30mm	40mm

## Open End



Ggauge	OD	length
27G	0.40mm	40mm



- \* We are able to handle anything from supply of our standard size injection needles on an OEM basis to customer requirements for specified dimensions.
- \* Manufacture of special sizes of injection needles is accommodated through sample manufacture.